

Rudolph's User Interface Software Version Compatibility with Windows 11

Rudolph Research has successfully installed and tested the following software versions on a PC using the Windows 11 operating system. All versions from that point on are considered compatible.

Autopol IV, V, V+, VI Polarimeters	3.2.6.1006 or later
Autopol I, II, III Polarimeters	2.2.5.0 or later
DDM Density Meters	3.9.9.0 or later
J457 Refractometers	4.5.2.1009 or later
Jx57 Refractometers	5.5.2.1009 or later
AutoFlex Automation	1.5.0.0 or later.

If additional information is required, please contact Rudolph Research Analytical.

Signature:



Name:

Bob Taggart

Title:

Vice President Engineering

Date:

16-Aug-2022